

Instantaneous 3D imaging of flame species using coded laser illumination

Elias Kristensson, Zheming Li, Edouard Berrocal, Mattias Richter and Marcus Aldén


Division of Combustion Physics, Lund University



LUND UNIVERSITY

Introduction


What is coded illumination?



Coded illumination


1 Superimpose structure

Guide light through a grid pattern




2 Modulated laser sheet

Sinus-modulated laser sheet through a cuvette



3 Camera view

View seen by a camera at 90 degrees




Coded illumination

2

Modulated laser sheet


Sinus-modulated laser sheet through a cuvette



3

Camera view


View seen by a camera at 90 degrees



4

Shift pattern


Distinguish between signal and background interferences



Coded illumination


3 Camera view

View seen by a camera
at 90 degrees




4 Shift pattern

Distinguish between signal
and background interferences



5 Final image


Calculate RMS on a
pixel level



Let's take a Fourier walk


1 Camera view

Sinus-modulated laser sheet through a cuvette




2 Fourier transform

Band-pass filter a small region (0.5% image power)



3 Filtered view

View the information stored locally in the FT



Need for high resolution?

1 2D signal




PLIF image of the OH radical

2 Zoomed view

Magnified view to see fine details

3 Fourier transform

Apply a low-pass filter on the image



B. Zhou, et al., Combustion and Flame (162) (2015) 2937-2953.

Simultaneous imaging

1

Signal 1

Simulated image of
a signal in 2D

2


Signal 2

Simulated image of a
signal in 2D


3

Both


Signal 1 + signal 2



+




=



Simultaneous imaging


2 Signal 2

Simulated image of a
signal in 2D




3 Both

Signal 1 + signal 2



4 Different
colors


Spectrally separate
signals



Simultaneous imaging


3 Both

Signal 1 + signal 2




4 Different colors

Spectrally separate signals



5 Time gating

Temporally separate signals




Simultaneous imaging

4

Different colors


Spectrally separate signals



5

Time gating


Temporally separate signals



6


Spatial modulation

Modulate each signal in space



Spatial lock-in detection


Theory



Spatial lock-in detection


1 Sample

Object structures that
we will probe




2 Modulation

An intensity modulated
illumination field



3 Modulated
signal

Modulated 2D image of
the sample



\times


$=$



Spatial lock-in detection


2 Modulation

An intensity modulated illumination field




3 Modulated signal

Modulated 2D image of the sample



4 Background


Not induced by the illumination



Spatial lock-in detection


3 Modulated signal

Modulated 2D image of the sample




4 Background

Not induced by the illumination



5 Without modulation


Detected signal without modulation



Spatial lock-in detection


4 Background

Not induced by the illumination




5 Without modulation

Detected signal without modulation



6 With modulation

Detected signal with modulation




Spatial lock-in detection

5

Without
modulation



Detected signal
without modulation



6

With
modulation


Detected signal with
modulation



7

Fourier
transform

Two peaks corresponding
to the modulation




Spatial lock-in detection

6

With
modulation


Detected signal with
modulation



7

Fourier
transform


Two peaks corresponding
to the modulation



8

Shift FT


Shift modulation peak to
origin



Spatial lock-in detection


7 Fourier transform

Two peaks corresponding to the modulation




8 Shift FT

Shift modulation peak to origin



9 Low-pass filter


Apply a low-pass filter on the shifted data



Spatial lock-in detection


8 Shift FT

Shift modulation peak to origin




9 Low-pass filter

Apply a low-pass filter on the shifted data



10 Inverse Fourier transform

Reveals the information carried by the modulation



Spatial lock-in for multiple detection

1

Signal 1

Sample 1 illuminated
with modulation

2


Signal 2

Sample 2 illuminated with
different modulation


3

Camera
view


Camera sees both
signals simultaneously



+




=



Spatial lock-in for multiple detection


2 Signal 2

Sample 2 illuminated with different modulation




3 Camera view

Camera sees both signals simultaneously



4 Fourier transform

Signal peaks are isolated from each other




Spatial lock-in for multiple detection

3

Camera view


Camera sees both signals simultaneously



4

Fourier transform


Signal peaks are isolated from each other



5

Inverse transform

Signal I extracted from inverse Fourier transform




Spatial lock-in for multiple detection

4

Fourier transform


Signal peaks are isolated from each other



5

Inverse transform


Signal 1 extracted from inverse Fourier transform



6


Inverse transform

Signal 2 extracted from inverse Fourier transform




Experimental setup


Modulated illumination for multiple exposures



Experimental setup




Camera view



First test


1 Cuvette

Cuvette with dye
illuminated with 2 sheets




2 Camera view

Metal piece in front of
one sheet



3 Top view


Sheets separated
~4 mm



First test


2 Camera view

Metal piece in front of one sheet




3 Top view

Sheets separated ~4 mm



4 Raw data


Image acquired with camera



First test


3 Top view

Sheets separated
~4 mm




4 Raw data

Image acquired with
camera



5 Magnified region


Modulation
frequencies different



First test


4 Raw data

Image acquired with camera




5 Magnified region

Modulation frequencies different



6 Fourier transform

Each signal as peaks in Fourier domain





First test

5

Magnified region



Modulation frequencies different



6

Fourier transform



Each signal as peaks in Fourier domain



7

Subimage 1

Information in band-pass filter #1




First test

6

Fourier transform


Each signal as peaks in Fourier domain



7

Subimage 1


Information in band-pass filter #1



8


Subimage 2

Information in band-pass filter #2



Results


Modulated PLIF images of formaldehyde distribution



Results


1 Camera view

4 superimposed images of CH₂O




2 Zoomed view

Different modulation patterns visible



3 Fourier transform

Signal components seen as peaks




Results

2

Zoomed view


Different modulation
patterns visible



3

Fourier transform


Signal components
seen as peaks



4

Decoded

Each signal
component decoded




Results

3

Fourier transform


Signal components
seen as peaks



4

Decoded


Each signal
component decoded



5

Distributed in 3D


Each signal is
positioned in space



Results


4 Decoded

Each signal
component decoded



5 Distributed in 3D

Each signal is
positioned in space



6 3D rendering


Instantaneous 3D
view of flame



4D imaging?


1 Droplets of dye

Dye taking the shape of a donut




2 “4D” view

Motion tracked in time



3 “4D” view

Motion tracked in time



Summary



Simultaneous detection
of several images

Images given different
structural “codes”

Spatial lock-in
demodulation




One camera solution

Perfect pixel-to-pixel
overlap

Potential for 4D imaging

Benefit from technological
development




Optical access

Trade resolution for
simultaneous detection

Signals share dynamic
range



Thank you for your attention!



Contact information:
email: elias.kristensson@forbrf.lth.se
tel: +46 46 222 4756

